SLE	REF N5402-REV 2	EN
Reusable Flow Sensor		
		T: +44 (0)208 681 1414 F: +44 (0)208 649 8570 E: <a href="mailto:info@sle.co.uk">info@sle.co.uk</a>

#### Instructions for Use

The sensor is intended for use with SLE4000, SLE5000 and SLE6000 ventilators. Refer to the ventilator user manual for ventilator specific connection and calibration instructions, cautions, warnings and specifications for the sensor.

The sensor is a reusable device. It is supplied clean but not sterile.

Before first use and after each subsequent use the flow sensor must be cleaned, disinfected or high level disinfected, and inspected as per the processing instructions. It can be reprocessed using the processing instructions up to 25 times or until it shows any signs of damage.

The flow sensor is used to detect the movement of gasses in the breathing circuit, which enables the ventilator to display and monitor the gas flow to and from the patient enabling monitoring of ventilation status.

Before insertion into the patient circuit the user must calibrate the flow sensor.

The flow sensor should be calibrated every 24 hours whilst in use, if the patient's condition permits.

**CAUTION:** The flow sensor may require cleaning during use.

**CAUTION:** Do not subject the sensor to high impact loads and shocks. This may destroy the sensor wires and render the flow sensor inoperable.

**WARNING:** Before use the sensor must be checked for damage. Damaged parts must not be used.

**WARNING:** Do not administer surfactant or use any nebulized gases (medications, salt solutions, etc.) in conjunction with the sensor, as they are likely to degrade the performance of the sensor and subsequent displayed accuracies.

#### Sensor Cleaning During Use

1-Disconnect the sensor body from the flow sensor cable.

2-Rinse the sensor body in sterile water with a gentle backwards and forwards motion.

3-Drain any water from the sensor.

4-Calibrate the flow sensor and re-insert sensor into patient circuit.

**WARNING:** Do not clean the flow sensor with compressed air or water jet. This may destroy the sensor wires and render the flow sensor inoperable.

#### Processing Instructions

**Cleaning:** Immediately after use rinse the sensor body in sterile water with a gentle backwards and forwards motion. A soap solution or mild alkaline solution may be used.

**Disinfection:** Use commercially available disinfectants that are recommended for use with PLASTIC MATERIALS. Immersion times and concentrations stated must be in accordance with manufacturer's instructions. Do not use disinfectants containing compounds similar to PHENOL or ALKYLAMINES (Glucorrotamine).


Eliminate all residues of cleaning agents and disinfectant used by thoroughly rinsing with sterile water after each cleaning and disinfecting procedure.

**High Level Disinfection:** The flow sensor is suitable for autoclaving. The sensor must not be connected to other standard connectors when autoclaved to prevent cracking. Ensure no other components/items are lying on the sensor during the autoclaving process.

134°C (277°F) (Allowable variation of temperature of +3°C) at 220kPa (32psi) with a minimum holding time of 3 minutes. or

121°C (248°F) (Allowable variation of temperature of +3°C) at 96kPa (14.1psi) with a minimum holding time of 15 minutes.

**Inspection:** Before each use the sensor must be checked for visible damage, cracking or soiling. Appropriate records should be maintained of the inspection results. Damaged sensors should be disposed of as clinical waste and replaced.

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