
**REF** L6000/SP2/001

## **SLE uSpO2 cable (Masimo SET)**

|           |                                  |
|-----------|----------------------------------|
| <b>EN</b> | <b>Instructions for use</b>      |
| <b>FR</b> | <b>Notice d'utilisation</b>      |
| <b>ES</b> | <b>Instrucciones de uso</b>      |
| <b>PT</b> | <b>Instruções de utilização</b>  |
| <b>IT</b> | <b>struzioni per l'uso</b>       |
| <b>NL</b> | <b>Gebruiksaanwijzing</b>        |
| <b>DE</b> | <b>Gebrauchsanleitung</b>        |
| <b>SE</b> | <b>Bruksanvisning</b>            |
| <b>EL</b> | <b>Οδηγίες χρήσης</b>            |
| <b>TR</b> | <b>Kullanım talimatları</b>      |
| <b>PL</b> | <b>Instrukcja użycia</b>         |
| <b>UA</b> | <b>Інструкція з використання</b> |
| <b>RU</b> | <b>Руководство пользователя</b>  |
| <b>JP</b> | <b>使用方法</b>                      |
| <b>CN</b> | <b>使用说明</b>                      |

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### SLE uSpO<sub>2</sub> Pulse Oximetry Cable (Masimo SET)

For use with Masimo sensors.

See sensor Instructions for use for sensor performance and use information.

#### Product description

The SLE uSpO<sub>2</sub> Pulse Oximetry Cable (Masimo SET) must be connected to the SLE6000 infant ventilator for power, communication, display, and alarm management.

#### Indications for use

Refer to the SLE6000 infant ventilator user manual for instructions for compatible sensors, accuracy information and cable connection.

The SLE uSpO<sub>2</sub> Pulse Oximetry Cable (Masimo SET) is indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor). The SLE uSpO<sub>2</sub> Pulse Oximetry Cable (Masimo SET) is indicated for use with neonatal, infant and paediatric patients up to 30kg during both motion and no motion conditions, and for patients who are well or poorly perfused in hospitals and hospital-type facilities.

#### Warnings

- As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.
- Do not place the Oximetry Cable or accessories in any position that might cause it to fall on the patient.
- Do not use the Oximetry Cable during magnetic resonance imaging (MRI) or in an MRI environment.
- Do not use the Oximetry Cable if it appears or is suspected to be damaged.
- Electrical Shock Hazard: to protect against injury, follow the directions below:
  - Avoid placing the device on surfaces with visible liquid spills.
  - Do not soak or immerse the device in liquids.
  - Use cleaning solutions only as instructed in this *Directions for Use*.
  - Do not attempt to clean the Oximetry Cable while monitoring patient.
- Do not place the device on electrical equipment that may affect the Oximetry Cable, preventing it from working properly.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- This device is only indicated for use with the SLE6000 and can only be connected to the SLE6000.
- The Oximetry Cable must be used with the SLE6000 infant ventilator.
- The Oximetry Cable is intended only as an adjunct device in patient assessment. It should not be used as the sole basis for diagnosis or therapy decisions. It must be used in conjunction with clinical signs and symptoms.
- When patients are undergoing Photo Dynamic therapy they may be sensitive to light sources. Pulse oximetry may be used only under careful clinical supervision for short time periods to minimize interference with Photo Dynamic therapy
- Displayed SpO<sub>2</sub> is empirically calibrated in healthy adult volunteers with normal levels of carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- Inaccurate SpO<sub>2</sub> readings may be caused by:
  - Elevated levels of COHb or MetHb: High levels of COHb or MetHb may occur with a seemingly normal SpO<sub>2</sub>. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
  - Intravascular dyes such as indocyanine green or methylene blue
  - Externally applied colouring and texture such as nail polish, acrylic nails, glitter, etc.
  - Elevated levels of bilirubin
  - Severe anemia
  - Low arterial perfusion
  - Motion artifact
- If SpO<sub>2</sub> values indicate hypoxemia, a laboratory blood sample should be taken to confirm the patient's condition.

## Serviceability


The SLE uSpO<sub>2</sub> Pulse Oximetry Cable has no serviceable parts. If damaged the cable should be replaced as it cannot be repaired.

## Service life

One year.


## Sensor Connection

Refer to Figure 1. Firmly insert the 9-pin sensor connector to the Oximetry Cable 9-pin receptacle shroud.

**Figure 1**

Refer to Figure 2. Rotate the clear latch over the mated connector until it snaps into place.

**Figure 2**

## Disconnection

To disconnect the sensor, reverse the above steps.

## Cleaning

The exterior surfaces of the Oximetry Cable may be cleaned with the following:

70% isopropyl alcohol

## Caution

Do not submerge the Oximetry Cable in any cleaning solution or attempt to sterilize by autoclave, irradiation, steam, gas, ethylene oxide or any other method. This will seriously damage the Oximetry Cable.

Disposal of Product: Comply with local laws in the disposal of the instrument and/or its accessories.

## Environmental

Operating Temperature

41°F to 104°F (5°C to 40°C)

Storage/Transport Temperature

-4°F to 122°F (-20°C to 50°C)

Operating Humidity

15% to 95% non-condensing

Storage/Transport Humidity

10% to 90% non-condensing

## Warranty

SLE warrants to the initial buyer only that each product it manufactures, when used in accordance with the directions provided with the Products by SLE, will be free of defects in materials and workmanship for a period of twelve (12) months.

The foregoing is the sole and exclusive warranty applicable to the products sold by SLE to buyer. SLE expressly disclaims all other oral, express or implied warranties, including without limitation any warranties of merchantability or fitness for a particular purpose. SLE's sole obligation and buyer's exclusive remedy for breach of any warranty shall be, at SLE's option, to repair or replace the product.

## Warranty exclusions

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

In no event shall SLE be liable to buyer or any other person for any incidental, indirect, special or consequential damages (including without limitation lost profits), even if advised of the possibility thereof. In no event shall SLE's liability arising from any products sold to buyer (under a contract, warranty, tort or other claim) exceed the amount paid by buyer for the lot of product(s) involved in such claim. In no event shall SLE be liable for any damages associated with a product that has been reprocessed, reconditioned or recycled. The limitations in this section shall not be deemed to preclude any liability that, under applicable products liability law, cannot legally be precluded by contract.

## No implied license




Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized sensors or cables which would, alone or in combination with this device, fall within the scope of one or more of the patents relating to this device.

## Caution:

For professional use. See sensor instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labelling.

| Symbol  | Definition                   |
|---|------------------------------|
|  | Caution                      |
|  | Consult Instructions for Use |
|  | Follow Instructions for Use  |

|  |   |
|--|---|
|    | Manufacturer  |
|   | Date of Manufacture   |
|   | Mark of Conformity to European Medical Device Directive 93/42/EEC   |
|   | Lot code  |
|   | Serial Number   |
|   | Catalogue number (model number)                                     |
|   | Storage Humidity Limitation   |
|   | Storage Temperature Limitation,                                     |
|   | Pressure Limitation   |
|   | Do not discard  |
|   | Separate collection for electrical and electronic equipment (WEEE). |
|   | Keep Dry  |
|  | Fragile, handle with care   |

Patents: <http://www.masimo.com/patents.htm>

For Compatibility Information Reference: [www.masimo.com](http://www.masimo.com)

Masimo, SET, are federally registered trademarks of Masimo Corporation.

## End-user license agreement

THIS SECTION of the DOCUMENT IS A LEGAL AGREEMENT BETWEEN YOU ("PURCHASER") AND SLE. IF YOU DO NOT AGREE TO THE TERMS OF THIS AGREEMENT, PROMPTLY RETURN THE ENTIRE PACKAGE, INCLUDING ALL ACCESSORIES, IN THEIR ORIGINAL PACKAGE, WITH YOUR SALES RECEIPT TO SLE FOR A FULL REFUND.

1. Grant of License. In consideration of payment of the license fee, which is part of the price paid for this product, SLE grants to Purchaser a non exclusive, nontransferable license, without right to sublicense, to use the copy of the incorporated software/firmware, and documentation in connection with Purchaser's use of the Masimo Products for their labeled purpose. SLE reserves all rights not expressly granted to Purchaser.

2. Ownership of Software/Firmware. Title to, ownership of, and all rights and interests in, any Masimo software and/or firmware and the documentation, and all copies thereof, remain at all times vested in Masimo Corporation, licensor to SLE, and they do not pass to Purchaser.

3. Assignment. Purchaser shall not assign or transfer this License, in whole or in part, by operation of law or otherwise, without SLE's prior written consent; any attempt without such consent, to assign any rights, duties or obligations arising hereunder shall be void.

4. Copy Restrictions. The software/firmware, mask works, circuit board layouts, and accompanying written materials are copyrighted. Unauthorized copying of the software, including software that has been modified, merged, or included with other software, or other written materials is expressly forbidden. You may be held legally responsible for any copyright infringement that is cause or incurred by your failure to abide by the terms of this license. Nothing in this license provides any rights beyond those provided by 17 U.S.C. §11 7.

5. Use Restriction. As the Purchaser, you may physically transfer the products from one location to another provided that the software/ firmware is not copied. You may not electronically transfer the software/firmware from the products to any other device. You may not disclose, publish, translate, release distribute copies of, modify, adapt, translate, reverse engineer, decompile, disassemble, or

create derivative works based on the Masimo Product, the software/firmware, or the written materials without the prior written consent of Masimo. Masimo Sensors that are designated for single use are licensed under Masimo patents for use on a single patient only, and are not sold. There is no license, implied or otherwise, that would allow use of single use Masimo Sensors beyond their intended single use. After use of single use Masimo Sensors, there is no further license granted by Masimo to use the sensors and they must be discarded.

6. Transfer Restrictions. The software/firmware is licensed to the Purchaser, and may not be transferred to anyone, except other end-users, without the prior written consent of SLE. In no event may you transfer, assign, rent, lease, sell, or otherwise dispose of the software/firmware or the products on a temporary basis.

7. Beneficiary. Masimo Corporation is a Beneficiary of this Agreement and has the right to enforce its provisions.

8. U.S. Government Rights: If you are acquiring software (including the related documentation) on behalf of any part of the United States Government, the following provisions apply: the software is deemed to be "commercial software" and "commercial computer software documentation," respectively pursuant to DFAR Section 227.7202 FAR I 2.212, as applicable. Any use, modification, reproduction, release, performance, display or disclosure of the software (including the related documentation) by the U.S. Government or any of its agencies shall be governed solely by the terms of this Agreement and shall be prohibited except to the extent expressly permitted by the terms of this agreement.

### Caution - Sample Pack

The sample pack supplied with this cable contains one sensor that is not to be used with the SLE uSpO2 cable (Masimo SET) and the SLE6000 infant ventilator as they are for patient weight sizes outside the intended use of the ventilator.

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All rights reserved. No part of this publication may be reproduced, stored in any retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopy, recording or otherwise, without prior permission of SLE. © Copyright SLE 11/07/2018
Tous droits réservés. Aucune partie de cette publication ne peut être reproduite, mise en mémoire dans un système de recherche documentaire ou transmise, sous quelque forme que ce soit ou par quelque procédé que ce soit, électronique, mécanique ou par photocopie, enregistrement ou autre, sans l'autorisation préalable de SLE. © Copyright SLE 11/07/2018.
Todos los derechos reservados. Queda prohibido reproducir o guardar esta publicación, en parte o en su totalidad, en cualquier medio de recuperación, así como comunicar su contenido de cualquier forma o por cualquier medio, ya sea electrónico, mecánico, fotocopia, grabación o de cualquier otra forma, sin el permiso previo de SLE. © Copyright SLE 11/07/2018.
Tutti i diritti riservati. Nessuna parte della presente pubblicazione può essere riprodotta, memorizzata in sistemi di recupero o trasmessa in qualsivoglia forma o con qualsiasi mezzo, elettronico, meccanico, mediante fotocopia, registrazione o altro, senza il preliminare consenso di SLE. © Copyright SLE 11/07/2018.
© Copyright SLE 11/07/2018
Alle Rechte vorbehalten. Ohne Erlaubnis von SLE darf kein Teil dieser Veröffentlichung in irgendeiner Form oder Weise (elektronisch, mechanisch, durch Fotokopieren, Aufnahme, etc.) reproduziert, auf Datenabfragesystemen gespeichert oder weitergegeben werden. © Copyright SLE 11/07/2018.