




SLE6000H, SLE6000C & SLE6000N

What is included with each new version?

Part Number	Mode of Ventilation	SLE6000H (Z6000H)	SLE6000C (Z6000C)	SLE6000N (Z6000N)
Invasive Ventilation				
Z6000/HFO	High Frequency Oscillation Ventilation (HFOV)	✔ <i>Option</i>	✘ <i>N/A</i>	✘ <i>N/A</i>
Standard	Conventional Mechanical Ventilation (CMV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Standard	Continuous Positive Airway Pressure (CPAP)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Standard	Synchronised Intermittent Mechanical Ventilation (SIMV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Standard	Synchronised Intermittent Mechanical with Pressure Support (SIMV+PS)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Standard	Pressure Supported Ventilation (PSV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Standard	Patient Triggered Ventilation (PTV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Z6000/VTV	Volume Targeting Ventilation (VTV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✘ <i>N/A</i>
Non-Invasive Ventilation				
Z6000/SLN	Nasal Continuous Positive Airway Pressure (nCPAP)	✔ <i>Standard</i>	✔ <i>Option</i>	✔ <i>Standard</i>
Z6000/NIP	Nasal Intermittent Positive Pressure Ventilation (NIPPV)	✔ <i>Standard</i>	✔ <i>Standard</i>	✔ <i>Standard</i>
Z6000/NIP	Triggered Nasal Intermittent Positive Pressure Ventilation (NIPPV Tr)	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Standard</i>
Standard with HFOV License	Nasal High Frequency Oscillation Ventilation (nHFOV)	✔ <i>Option</i>	✘ <i>N/A</i>	✘ <i>N/A</i>
Z6000/SLN	DuoPap	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Standard</i>
Z6000/O2T	O ₂ Therapy	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Standard</i>
Z6000/CLP	OxyGenie® License	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Standard</i>
Z6000/SPO and L6000/SP2	SpO ₂ Sensor and License	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Option</i>
Z6000/ETC and LETC2/KIT	CO ₂ Module and License	✔ <i>Option</i>	✔ <i>Option</i>	✘ <i>N/A</i>
Z6000/ITB	IntelliBridge	✔ <i>Option</i>	✔ <i>Option</i>	✔ <i>Option</i>
Standard	RS232 (v2.0 and v3.0)	✔ <i>Standard</i>	✔ <i>Standard</i>	✔ <i>Standard</i>



Standard in selected ventilator model Z6000H Z6000C Z6000N.

Option – accessory and license to be selected.

N/A not available in selected main unit.