

First Breath™ Bonnet

Access to scalp

Attachment block

Single Tube


Ribbon ties

Nasal Prong or Mask

Hook-and-Loop tab

- Available in 12 colour-coded sizes
- Hook-and-Loop tab designed to allow for adjustment
- Enables access to the scalp
- Single point of fixation designed to aid alignment
- Attachment block designed to transfer weight from the nasal area
- Compatible with First Breath nCPAP Generator and has two fixation options - the First Breath Headgear and the new First Breath Bonnet

First Breath Headgear



First Breath™ Bonnet

First Breath Bonnet Ordering Information

First Breath Bonnet - White (16-18cm)	VXWH000/10
First Breath Bonnet - Grey (18-20cm)	VXGY00/10
First Breath Bonnet - Pink (20-22cm)	VXPK0/10
First Breath Bonnet - Brown (22-24cm)	VXBRI/10
First Breath Bonnet - Yellow (24-26cm)	VXYE2/10
First Breath Bonnet - Blue (26-28cm)	VXBL3/10
First Breath Bonnet - Gold (28-30cm)	VXGO4/10
First Breath Bonnet - Green (30-32cm)	VXGR5/10
First Breath Bonnet - Burgundy (32-34cm)	VXBU6/10
First Breath Bonnet - Orange (34-36cm)	VXOR7/10
First Breath Bonnet - Dark Green (36-38cm)	VXDG8/10
First Breath Bonnet - Navy (38-40cm)	VXNA9/10

MK-PHL-000006

This document has been produced by Inspiration Healthcare Group.



Manufactured by

Viomedex Limited, Unit 12 - 14 Swan Barn Business Centre Old Swan Lane Hailsham East Sussex BN27 2BY

+44 (0)330 175 0000 | info@inspiration-healthcare.com | inspirationhealthcaregroup.com

The details given in this brochure are correct at time of going to press. Any specifications are subject to change and may be different depending on where in the world you are located. This information does not supersede or replace the Instructions for Use. Always refer to the latest version of Instructions For Use prior to using any products.