



InspireView™

Video Laryngoscope for Preterm and Term Infants

InspireView™

In the delicate and precise field of neonatal care, every moment is critical. InspireView™ transforms the way clinicians visualise and manage the airway. This can help in allowing endotracheal tubes to be placed correctly and while providing a clear view of the airway structures, making intubation safer and more efficient.

InspireView's compact, three-part design allows for rapid setup, eliminating the need for cumbersome cables and ensuring mobility. From its suitability for the tiniest infants to its easy-to-clean reusable parts and single-use blades, InspireView is engineered to enhance uptime and streamline your workflow.

It's as easy as 1, 2, 3...

1

Touch Screen

- The 3" touchscreen provides a clear image
- Auto shut down after 10 mins of inactivity
- 3 hour rechargeable battery
- Fast access to photos and videos

2

Reusable Handle

- Allows you to capture clear photos and videos with thumb controlled button

3

Disposable Design

- Polycarbonate single patient use blade and handle.
- A selection of blades for preterm and term infants, allowing a clear view of the glottis



Other Benefits

Visualise the patient airway naturally and effectively.

Illumination and high-resolution views of airway structures.

Set up quickly with our compact, three-part device, perfect for urgent situations.

Experience the freedom of a cable-free design, enhancing mobility and convenience.

Utilise blades and handles specifically designed for the smallest infants.

Benefit from easy and fast cleaning of reusable parts and single-use blades, ensuring more time for patient care.

Capture videos and photos for training purposes or patient records, enhancing educational opportunities and documentation.

Ordering Information

InspireView Video Laryngoscope System comprising of InspireView Monitor with Touchscreen and Handle	UET-A0
InspireView Video Laryngoscope Monitor and Touchscreen for use with InspireView Video Laryngoscope Handle	UED-M3-S
InspireView Video Laryngoscope Handle for use with InspireView Video Laryngoscope monitor and Touchscreen	UET-A0 Handle
InspireView Single Use Miller Blade Unit for Preterm Infants (Size 0)	UED-A0
InspireView Single Use Miller Blade Unit for Term Infants (Size 1)	UED-MA1
InspireView Single Use Curved Blade Unit for Term Infants (Size 1)	UED-A1
InspireView Carry Case	UED-Case

Manufactured by

Zhejiang UE Medical Corp
8 Youyi Rd, Baita Economic Development Zone,
Xianju, Zhejiang 317300,
China
T +86-571-28057363
E export@ueworld.com
W ueworld.com

Distributed by



Commerce Park,
Commerce Way
Croydon, CR0 4YL
England
T +44 (0)330 175 0000
E info@inspiration-healthcare.com
W inspirationhealthcaregroup.com