


SLE ESSENTIAL


High Flow Nasal Cannula




Where Simplicity Meets Purpose

For clinicians seeking High Flow on their SLE Ventilators, the new SLE High Flow Nasal Cannula can help deliver precise flow for an **affordable price**.

Designed for **simplicity** and **flexibility**, it supports efficient workflows while adapting to individual patient needs.



Features & Benefits



Affordable Access

Helps manage costs for budget-conscious facilities



Simple Setup

Streamlined design for ease of use



Flexible Fit

Adaptable fixation options for patient comfort

Technical Specifications & Ordering Information

Nasal Cannulas

	 VINC-01	 VINC-02	 VINC-03	 VINC-04
Size	S	M	L	XL
Part Number	VINC-01	VINC-02	VINC-03	VINC-04
Lock Ring Colour	Uncoloured	Green	Blue	Orange
Nasal Prongs OD	2.8 mm	2.8 mm	3.7 mm	4.5 mm
Space Between the Two Nasal Prongs	2.5 mm	4.2 mm	4.8 mm	5.1 mm
Max. Flow	15 l/min	15 l/min	25 l/min	25 l/min

High Flow Nasal Circuits

Part Number	Description
C4330	High Flow Nasal Cannula Circuit
C4331	High Flow Nasal Cannula Circuit with Autofeed Chamber

MK-PHL-000116



Distributed by

SLE Limited

Commerce Park, Commerce Way,
Croydon, CR0 4YL

England

T +44 (0)330 175 0000

E info@inspiration-healthcare.com

W inspirationhealthcaregroup.com